

## Activity 3.4

### Investigate and describe the effects on plant tissue of immersing them in different solutions

#### Skills

A03.1 Using techniques, apparatus and materials

A03.3 Observing, measuring and recording

A03.4 Interpreting and evaluating observations and data

#### Safety

Take care when using a sharp blade to cut the plant tissue.

- 1 Set up a microscope.
- 2 Take three clean microscope slides. Label them A, B and C.
- 3 Put a drop of distilled water onto the centre of slide A.
- 4 Put a drop of medium concentration sugar solution onto slide B.
- 5 Put a drop of concentrated sugar solution onto slide C.
- 6 Peel off a very thin layer of coloured epidermis from a *Rhoeo* leaf, or other leaf, or from a rhubarb petiole. To get good results, it should be as thin as possible (only one cell thick).
- 7 Cut three squares of this epidermis, each with sides about 5 mm long,
- 8 Put one square into the drop of solution on each of your three slides.
- 9 Carefully cover each one with a coverslip. Clean excess liquid from your slides with filter paper.
- 10 Look at each of your slides under the microscope. Make a labelled drawing of a few cells from each one.



#### Questions

- A1 Which part of the cell is coloured?
- A2 What has happened to the cells in pure water? Explain your answer.
- A3 What has happened to the cells in medium concentration sugar solution? Explain your answer.
- A4 What has happened to the cells in concentrated sugar solution? Explain your answer.